



قسم الكيمياء الحيوية  
Biochemistry Department  
College of Science - King Saud University

## **BCH 471**

### Experiment (5)

# **Determination of plasma enzymes**

Determination of LDH in serum

# OBJECTIVES

- To determine the level of LDH in serum.
- To evaluate the presence of tissue damage.

# PLASMA ENZYMES

- Plasma enzymes are classified into:
  1. Functional plasma enzymes
  2. Non functional plasma enzyme

## *Differences between functional and non-functional plasma enzymes*

	<b>Functional plasma enzymes</b>	<b>Non-functional plasma enzymes</b>
<b>Concentration in plasma</b>	Present in plasma in higher concentrations in comparison to tissues	Normally, present in plasma in very low concentrations in comparison to tissues
<b>Function</b>	Have known functions	No known functions
<b>The substrates</b>	Their substrates are always present in the blood	Their substrates are absent from the blood
<b>Site of synthesis</b>	Liver	Different organs e.g. liver, heart, brain and skeletal muscles
<b>Effect of diseases</b>	Decrease in liver diseases	Different enzymes increase in different organ diseases
<b>Examples</b>	Clotting factors e.g. prothrombin, Lipoprotein lipase and pseudo-choline esterase	ALT, AST, CK, LDH, alkaline phosphatase, acid phosphatase and amylase,

# SOURCES OF NON FUNCTIONAL PLASMA ENZYME

- Cell damage with the release of its content of enzymes into blood e.g. Myocardial infarction and viral hepatitis
- Obstruction of normal pathways e.g. Obstruction of bile duct increases alkaline phosphatase
- Increase of the enzyme synthesis e.g. bilirubin increases the rate of synthesis of alkaline phosphatase in obstructive liver disease
- Increased permeability of cell membrane as in hypoxia

# MEDICAL IMPORTANCE OF NON FUNCTIONAL ENZYMES


- Measurement of non functional enzymes is important for:
  - Diagnosis of diseases
  - Prognosis of the disease: following up of the treatment by measuring plasma enzymes before and after treatment.

# LACTATE DEHYDROGENASE (LDH)

- Lactic acid dehydrogenase (LDH) is an enzyme that helps produce energy.
- LDH is most often measured to evaluate the presence of tissue damage (diagnostic)
- The enzyme LDH is in many body tissues, especially the heart, liver, kidney, skeletal muscle, brain, blood cells, and lungs.

# LDH REACTION

- LDH is a hydrogen transfer enzyme which catalyzes the interconversion of pyruvate and lactate.
- Lactate is released into the blood and is eventually taken up by the liver.
- In the liver, it catalyzes the oxidation of L-lactate to pyruvate (L→P) with the mediation of NAD as hydrogen acceptor, eventually converting pyruvate to glucose.
- Glucose is released into the blood and then taken up by resting muscles, red blood cells, and other tissue.





- Exercising muscles convert (and red blood cells metabolize) glucose to lactate.
- The optimum pH for lactate pyruvate ( $L \rightarrow P$ ) reaction is 8.8 – 9.8
- While for pyruvate to lactate ( $P \rightarrow L$ ) is 7.7 – 7.8.
- The enzyme is inhibited by sulfhydryl reagents such as P-chloromercuribenzoate and mercuric ions.

# LDH ISOENZYMES

- LDH exists in 5 forms (isoenzymes), which differ slightly in structure.
- All of these isoenzymes can be measured in the blood, and can be separated by electrophoresis.

LDH Isoenzyme	Tissues
LDH-1	is found primarily in heart muscle and red blood cells.
LDH-2	is concentrated in white blood cells.
LDH-3	is highest in the lung
LDH-4	is highest in the kidney, placenta, and pancreas
LDH-5	is highest in the liver and in skeletal muscle


**↑ LDH  
in  
plasma**

Diseases	Examples
Myocardial infarction	
Liver Disease	Toxic jaundice
	Viral hepatitis
	Obstructive jaundice
Anemia	Pernicious anemia
	Megaloblastic anemia
Renal Diseases	Tubular necrosis
	Pyelonephritis
Malignant Disease	Lung Cancer
	Hodgkin's disease

# -Lactate Dehydrogenase Assay

## Principle:

LDH catalysis the following reaction:



The rate of NADH formation is indicated by **increase the absorbance** at 340nm and it is **directly proportional to serum LDH activity**.

If:

NADH is **product** : **increase** the absorbance /min

NADH is **reactant**: **decrease** the absorbance /min

# Method

	<b>Tube</b>
LDH reagent	3ml
Pre- warm at 37 c for 3 minutes	
Sample (serum)	0.1 ml(100 $\mu$ l)
Mix and incubate at 37 c for 1 min . Read the absorbance at 340 nm against distilled water every minute for 3 minute Determine $\Delta A/\text{min}$	

# RESULTS

	Time (min)	Absorbance at 340 nm
A1	1	
A2	2	
A3	3	

# CALCULATIONS

$$\Delta A_1 = A_2 - A_1$$

$$\rightarrow \Delta A/\text{min} = (\Delta A_1 + \Delta A_2) / 2$$

$$\Delta A_2 = A_3 - A_2$$

$$\text{LDH (U/L)} = \Delta A \times 4984$$

- **Normal Values** 109 to 245 U/L