

ESTIMATION OF INORGANIC PHOSPHATE IN MILK AND SOFT DRINK



OBJECTIVE

- ▶ Estimation of organic phosphate in milk and soft drink


PHOSPHATE IN FOOD

- ▶ Phosphate occurs naturally in the form of organic esters in many kinds of food, including meat, potatoes, bread, and milk.
- ▶ Phosphate also used as a food additive (inorganic phosphate) as a preservative, a flavor or color enhancer, extend shelf life, and retain moisture..

- ▶ Soft drinks are complex mixtures containing a variety of substances such as colouring compounds, flavoring agents, acidifiers, sweeteners, preservatives, and caffeine.
- ▶ The most common acidifier used in soft drinks is phosphoric which gives a tangy taste in the mouth.
- ▶ Phosphoric acid can also acts as a preservative, keeping the contents of the bottle fresh.



- Due to the use of phosphoric acid, cola is actually more acid than vinegar which no body can drink straight. But a ton of sugar, dyes and flavoring are added to mask the acidity.



PRINCIPLE:

- ▶ Phosphoric acid is colorless, they cannot be directly determined using visible-light spectrophotometry
 - ▶ Instead, we will quantitatively convert them into a colored substance, whose absorbance can be easily measured
1. Inorganic phosphate reacts with ammonium molybdate in an acid solution to form phosphomolybdic acid
 2. phosphomolybdic acid is then reduced by a reducing agent (3% ascorbic acid) to give molybdenum blue a green/ blue color but does not affect the uncombined molybdic acid .

METHOD

	Standard	Milk sample	Soft drink sample	Ammonium molybdate	Ascorbic acid
Blank	----	---	----	0.5 ml	0.5 ml
S1	2	---	----		
S2	2	---	-----		
S3	2	---	---		
S4	2	---	---		
S5	2	---	---		
M1	---	2	---		
M2	---	2	---		
SD1	---	--	2		
SD2	---	--	2		

METHOD

- ▶ Mix thoroughly after each addition .
- ▶ Allow to stand for 10 min
- ▶ (a deep blue/green colour should develop).
- ▶ Measure the absorbance at 650 nm.

- RESULTS AND CALCULATIONS:

- ▶ Plot a graph between absorbance and concentration of phosphate in various standard solutions and obtain the standard curve.
- ▶ From the curve determine the amount of phosphate in the test solution.