



قسم الكيمياء الحيوية  
Biochemistry Department  
College of Science - King Saud University

**BCH 471**

Experiment (4)

**Haemolysing Agents &  
Detection of blood**

# OBJECTIVES

- ✓ To find the nature of substances that cause haemolysis of RBCs.
- ✓ To detect the presence of blood.

# HAEMOLYSIS


- It is the breaking open of RBCs and the release of hemoglobin and the red cell contents into the surrounding fluid (plasma).
- The concentration of potassium inside red blood cells is much higher than in the plasma  
→ so elevated potassium is usually found in biochemistry tests of hemolysed blood.
- Haemolysis may occur in:
- Vivo:
  - *Causes:* Immune reactions, Infections, Medications, Toxins and poisons, or by a defect in the RBC membrane.
- Vitro:
  - *Causes:* Placing RBCs in a hypotonic solution.

# IMPORTANCE OF HAEMOLYSIS PRACTICALLY

- Breaking down RBCs to release their content is often necessary in biochemistry for:
  - Estimation of haemoglobin.
  - To obtain erythrocytes-free preparation of leukocytes and platelets.

# OSMOTIC PRESSURE


- Osmosis : It is the diffusion of water across a selectively permeable membrane into a region of higher solute conc., once an equilibrium is reached the flow of water stops .
- Osmotic pressure : It is the pressure that generates from osmosis .



(a) In osmosis, there is net movement of solvent from the area of **higher solvent concentration** (*lower solute concentration*) to the area of **lower solvent concentration** (*higher solute concentration*).


# ISOTONIC SOLUTION

- It is a solution that has the same salt concentration as the normal cells of the body and the blood, having equal osmotic pressure.
- For example: 0.9% NaCl solution, have the same osmotic pressure as serum → it does not affect the membranes of the RBCs.
- In hospitals, intravenous fluids are isotonic.




# HYPOTONIC SOLUTION

- In a hypotonic solution, there is a lower concentration of solute outside a cell, creating an environment with lower osmotic pressure than what is contained within the cell.
- As a result, the RBCs would hemolyze or burst.
- Examples of Hypertonic Solutions
  - 0.45% NaCl ( Lower concentration than normal saline 0.9%NaCl )
  - Distilled water
  - Dextrose 2 or 2.5% in water.



# HYPERTONIC SOLUTION

- In a hypertonic solution, there is a higher concentration of solute outside a cell, creating an environment with higher osmotic pressure than what is contained within the cell.
- As a result, the RBCs would shrink.
- Example of Hypertonic Solutions
  - 1.2% NaCl ( higher concentration than normal saline 0.9%NaCl )





# HOW TO CALCULATE THE OSMOLARITY OF A SOLUTION

**For example: Calculate the osmolarity of NaCl ?**

Osmolality of RBC = 0.308 Osmolar

$O = M \times \text{no. of dissociation particles}$

$\text{NaCl} \rightarrow \text{Na}^+ + \text{Cl}^-$  (no. of dissociation particles = 2)

$M = \text{no. of moles / volume (L)}$

$0.308 = M \times 2 \rightarrow M = 0.154 \text{ M}$

\*\*To calculate in w/w% expression

$M = \text{no. of moles} / V \text{ (in L)}$

no. of moles = weight / Molecular weight So,

weight =  $M \times \text{Molecular weight}$  (in 1000 ml)

$= 0.154 \times 58.5 = 9 \text{ g}$

$9 \text{ g} \rightarrow 1000 \text{ ml}$

$? \text{ g} \rightarrow 100 \text{ ml}$

w/w% of NaCl =  $(9 \times 100) / 1000 = 0.9\%$   $\rightarrow$  isotonic

# METHOD


- Into 6 dry clean test tubes (A, B, C, D, E, F), pipette 3 drops of the suspended RBC's in saline solution, and add to each tube as indicated the following table:

	A	B	C	D	E	F
NaCl 0.45%	5ml					
NaCl 1.2%		5ml				
Sucrose 6%			5ml			
NaOH 0.1M				3drops		
HCl 0.1M					3drops	
Dis. H <sub>2</sub> O						5ml
NaCl 0.9%				5ml	5ml	

- Wait 30 minutes.
- Observe whether Haemolysis has taken place

Mwt: NaCl= 58.5      Sucrose= 342.3

# RESULTS



Note that the hemolyzed sample is transparent, because there are no cells to scatter light.

# DETECTION OF BLOOD BY BENZIDINE TEST

- Principle:
  - This method depend on the fact that the haem group of haemoglobin possesses a peroxidase-like activity which catalyses the breakdown of hydrogen peroxide ( $H_2O_2$ ).
  - The oxidising species formed in this reaction can then react with benzidine giving a blue color.
- It is often necessary to detect the presence of small quantities of blood in urine, stomach contents etc.
- **The test is not specific for blood** because peroxidases are present in milk, potatoes and pus, as well as the ions of  $Fe^{+3}$ ,  $Cu^{+2}$  and  $K^+$  will give **false positive results..**

# METHOD

- 3ml of suspended blood cells solution is boiled in water bath for 3 minutes and then cool it under tap water.
- Add 2 ml of benzidine solution, followed by 1 ml of hydrogen peroxide solution.

# NOTE

The general principle is that:

- if the test is negative, blood is absent.

But

- if the test is positive, blood is probably, not definitely present.

→ For this reason the tests are often described as  
*"presumptive tests"*.