




Qualitative tests of Carbohydrate



Carbohydrate :

- Carbohydrates are the **key source of energy** used by living things.
- Also serve as extracellular structural elements as in cell wall of bacteria and plant.
- Carbohydrates are defined as the **polyhydroxy aldehydes or polyhydroxy ketones.**
- Most , but not all carbohydrate have a formula $(CH_2O)_n$ (hence the name hydrate of carbon)
- In human body, the D-glucose is used.
- Simple sugars ends with -ose.



- Several **classifications of carbohydrates** have proven useful, and are outlined in the following table.

| | | | | |
|---------------------|--|--|--|--|
| Complexity | Simple Carbohydrates monosaccharides | Complex Carbohydrates disaccharides, oligosaccharides & polysaccharides | | |
| Size | Tetrose C ₄ sugars | Pentose C ₅ sugars | Hexose C ₆ sugars | Heptose C ₇ sugars etc. |
| C=O Function | Aldose sugars having an aldehyde function or an acetal equivalent. Ketose sugars having a ketone function or an acetal equivalent. | | | |
| Reactivity | Reducing sugars oxidized by Tollens' reagent (or Benedict's or Fehling's reagents). Non-reducing sugars not oxidized by Tollens' or other reagents. | | | |

Classification :



1-Simple sugar (one unit) :

Monosaccharides contain **one** monosaccharide unit.

2-Complex sugar (more than one) :

- **Disaccharides** contain **two** monosaccharide units.
- **Oligosaccharides** contain **3-9** monosaccharide units.
- **Polysaccharides** can contain more than 9 monosaccharide units.


-Complex carbohydrates can be broken down into smaller sugar units through a process known as **hydrolysis**.




Amylose

Reducing and non reducing sugars

- Reducing and non reducing sugar :If the oxygen on the anomeric carbon of a sugar is **not attached** to any other structure, that sugar can act as a reducing agent and is termed a reducing sugar.



reducing



Non-reducing

Solubility of sugars [physical property]:

- Monosaccharide and disaccharide **can** be dissolved freely in water because water is a polar substance, while polysaccharide **cannot** be dissolved easily in water, because, it has **high molecular weight** , which give colloidal solutions in water.



Chemical Properties of Carbohydrates:

- 1- **Molisch Test**: specific for carbohydrates.
- 2- **Benedict's Test**: presence of reducing sugars.
- 3- **Barfoed's Test**: test used for detecting the presence of monosaccharides.
- 4- **Bial's Test**: used to detect pentose [5C] monosaccharides.
- 5- **Seliwanoff's Test**: distinguish between aldoses and ketoses.

1. Molisch test:


This test is specific for all carbohydrates Monosaccharide gives a **rapid** positive test, Disaccharides and polysaccharides react **slower**.

Objective:

To identify the carbohydrate from other macromolecules, lipids and proteins.


Principle:

- The test reagent(H₂SO₄) dehydrates pentose to form **furfural** and dehydrates hexoses to form **5- hydroxymethyl furfural**.
- The furfural and 5- hydroxymethyl furfural further react with α-naphthol present in the test reagent to produce a **purple ring**.



$\xrightarrow{\alpha\text{-naphthol}}$
 [Present in the reagent]

Purple ring



$\xrightarrow{\alpha\text{-naphthol}}$

Purple ring

Method :

- 1-Two ml of a sample solution is placed in a test tube.
- 2-Two drops of the Molisch reagent (which α -naphthol in 95% ethanol) is added.
- 3-The solution is then poured slowly into a tube containing two ml of concentrated sulfuric acid so that two layers form, producing violet ring appear as liaison between the surface separations.

| Tube | Observation |
|---------|-------------|
| Glucose | |
| Lactose | |
| Starch | |



2. Benedict's test:


Objective:

To detect the presence of reducing sugars.

- **All monosaccharides** are reducing sugars; they all have a free reactive carbonyl group.
 - **Some disaccharides** have exposed carbonyl groups and are also reducing sugars. Other disaccharides such as sucrose are non-reducing sugars and will not react with Benedict's solution
- Large polymers** of glucose, such as starch, are **not** reducing sugars, since the concentration of hemiacetal groups is very low.

Principle :


- The copper sulfate (CuSO₄) present in Benedict's solution reacts with electrons from the aldehyde or ketone group of the reducing sugar in **alkaline medium**.
- Reducing sugars are **oxidized** by the copper ion in solution to form a carboxylic acid and a **reddish precipitate** of copper (I) oxide.



Reducing sugar

carboxylic acid

reddish
precipitate
[copper oxide].



Method :

- One ml of a sample solution is placed in a test tube.
- Two ml of Benedict's reagent is added.
- The solution is then heated in a boiling water bath for five minutes.
- A positive test is indicated by: The formation of a reddish precipitate.

| Tube | observation |
|-----------|-------------|
| 1-glucose | |
| 2-lactose | |
| 3-starch | |

3.Barfoed's Test:

Objective :


To distinguish between mono- , di- and poly saccharides.

Principle :

Barfoed's test used copper (II) ions in a **slightly acidic medium**.

Reducing **monosaccharides** are oxidized by the copper ion in solution to form a carboxylic acid and a **reddish** precipitate of copper (I) oxide within three minutes. Reducing **disaccharides** undergo the same reaction, but do so at a **slower** rate.

-The **nonreducing** sugars give negative result.



Method :

- Place one ml of a sample solution in a test tube.
- Add 3 ml of Barfoed's reagent (a solution of cupric acetate and acetic acid).
- Heat the solution in a boiling water bath for 6 minutes(after the 3 min check the tubes).



| Tube | observation |
|---------|-------------|
| glucose | |
| lactose | |
| starch | |



4. Bial's Test:

Objective: To distinguish between pentose monosaccharide and hexose monosaccharide

Principle: Bial's test uses concentrated HCl as a dehydrating acid and orcinol + traces of ferric chloride [FeCl₃] as condensation reagent. The test reagent dehydrates **pentoses** to form furfural. Furfural further reacts with orcinol and the iron ion present in the test reagent to produce a **bluish or green** product, while **hexoses** yield muddy-brown to grey condensation product.



Method :

- Put 2 ml of a sample solution in a test tube.
- Add 3 ml of Bial's reagent (a solution of orcinol, HCl and ferric chloride) to each tube.
- Heat the tubes gently in hot water bath for two minutes.
- If the color is not obvious, more water can be added to the tube.

| Tube | observation |
|-----------|-------------|
| 1-glucose | |
| 2-ribose | |

5.Seliwanoff's Test:


Objective:

used to distinguish between aldoses (like glucose) and ketoses (like fructose).

Principle:

Seliwanoff's Test uses 6M HCl as dehydrating agent and resorcinol as condensation reagent. The test reagent dehydrates **ketohexoses** to form 5-hydroxymethylfurfural. 5-hydroxymethylfurfural further condenses with resorcinol present in the test reagent to produce a **cherry red product** within two minutes.

-**Aldohexoses** react to form the same product, but do so more slowly giving yellow color.



Method :

- One half ml of a sample solution is placed in a test tube.
- Two ml of Seliwanoff's reagent (a solution of resorcinol and HCl) is added.
- The solution is then heated in a boiling water bath for two minutes.

| Tube | observation |
|------------|-------------|
| 1-glucose | |
| 2-fructose | |

Questions:

- 1- Name the complex formed by the addition of concentrated sulfuric acid to sugar solution and explain the reaction?
- 2- Why sucrose gives negative Benedict test?
- 3- Explain, although starch has free hemiacetal bond it gives negative Benedict test?
- 4- Why glucose (monosaccharide) and lactose (disaccharide) give positive Benedict test?
- 5- What is the difference between Benedict and Barfoed's reaction?
- 6- What are the carbohydrates' that give positive result with Seliwanoff ? why?