

# Harvest and Staining Protocols


By: Sahar AlSubaie

# Harvest Protocols:

- ❑ Arresting cells at metaphase stage .
- ❑ Hypotonic treatment.
- ❑ Cells fixation.

# metaphase stage:

chromosomes are lined at the equatorial plate before sister chromatids are pulled to opposite poles by spindle fibers for incorporation into the two daughter cells.



# Arresting cells at metaphase stage .

Done by adding mitotic arrestant such as **Colcemid**.

# Mechanism of Colcemid action:

preventing spindle fibers formation and thereby preventing pulling the two sister chromatides to the daughter cells.

# Hypotonic treatment:

Using hypotonic solution results in :

- ✓ Swelling of the cells :  
Thinning cytoplasm & enhance chromosome spreading
- ✓ Lysis of RBCs.

# Cells fixation:

- ❑ Fixative is added to kill and preserve the cells.
- ❑ It removes the water from the cell (dehydrate them).
- ❑ hardens membranes and chromatin.

# Cells fixation:

- ❑ 1<sup>st</sup> few drops of fixative stop the action of hypotonic solution and prepare cells for next higher fixative concentrations
- ❑ The brown supernatant is the color of the methemoglobin. (formed after the addition of fixative)
- ❑ Multiple fixative steps should be carried out. (until a clear supernatant is obtained)



# Procedure:

- ❑ after incubation , add 50 µl of colcemid.
- ❑ Incubate for 15 mins at 37 °C .
- ❑ Centrifuge the tube for 8 mins at 1200 rpm speed .
- ❑ Discard the supernatant without touching the buffy coat ( because cells are there )
- ❑ resuspend pellets :Mix thoroughly using the your palm hand , continue until the buffy coat disappears.
- ❑ Add 2 ml of pre-warmed hypotonic solution(0.075 M KCl) drop by drop while mixing using vortex .
- ❑ Add 8 ml of hypotonic solution without mixing .
- ❑ Incubate at 37 °C for 15 mins.
- ❑ 1<sup>st</sup> wash : Add 1 ml of fixative(3:1 methanol: acetic acid) dropwise
- ❑ Centrifuge for 8 min at 1200 rpm.

# Procedure:

- Remove the supernatant but not completely .
- Mix thoroughly as the cells are so sticky
- Add 2 ml of fixative drop by drop using vortex .
- Add more fixative ( up to 10 ml) without using vortex.
- Centrifuge for 8 min at 1200 rpm.
- Now the 2<sup>nd</sup> wash: starts after removing the supernatant.
- Mix to detach the cells from the tube wall
- Add 10 ml of fixative
- Centrifuge for 8 min at 1200 rpm.
- Don't remove the supernatant now.
- Keep at 4 °C until you can do the dropping and staining

# Dropping Protocol:

# Procedure:

1. Replace the fixative with a new one. Using a fine glass pasture pipette make a gentle mix without making air bubbles, as this will cause cells to clump.
  2. Clean the slides with methanol, that will help to get a uniformity in spreading.
  3. Carry the slide in 45 degrees and drop the sample from about one meter apart.
  4. As a test slide drop about 4 drops on the slide, leave to air dry completely and observe the mitotic index and metaphases compactness.
- Aging: Put the slides in the oven at 60C over night or at 90C for 90 minutes, that will driving off water and get better banding pattern.

Staining :

# Staining Protocols:

- ❑ G banding (Giemsa)
- ❑ Q banding (Quinacrine, fluorescent stain )
- ❑ R banding ( Reverse)
- ❑ C banding (Centromeric(heterochromatin))
- ❑ Ag-Nor stain (Nucleolar Organizing Regions)

# G banding :

- ❑ The most common .
- ❑ Trypsin hydrolyses the protein component of the chromatin (Wang and Federoff) :  
allowing the Giemsa dye to react with the exposed DNA.

# G banding :

Giemsa stain is a complex mixture of dyes:

**The basic dye:** aminophenothiazin dyes azure A, azure B, azure C, thionin, and methylene blue.

**The acidic dye:** eosin.



# G banding :

Immersing slides in fetal calf serum to stop the activity of the trypsin since the serum contains  $\alpha$ 1-antitrypsin

# G banding :

- ❑ Appropriately stained chromosomes  
neither too dark nor too pale.
- ❑ Under-trypsinized chromosomes  
have indistinct bands and little contrast. i.e.fuzzy in appearance.

## G banding :

- ❑ Over-trypsinized chromosomes

have sharp bands and often appear frazzled at the ends.

- ❑ Eventually overtrypsinized chromosomes

are very pale after staining and may appear ghost-like and very swollen.

# G banding :

- ❑ Each chromosomes characterized by dark & light bands
- ❑ 400 bands in the haploid genome
- ❑ Light bands are genetically active sites
- ❑ Dark bands are gene poor rich in AT pairs.

## G banding :

- ❑ Chromosomes 13,18,21 gene poor (very dark chromosomes)
- ❑ Chromosome 21 is smaller than Chromosome 22
- ❑ Chromosome 21 (200 genes) is half as many as Chromosome 22 (400 genes)


# Procedure :

- Slides must cool down to RT
- Immerse in trypsin solution 10-15 sec
- Wash in phosphate buffer or serum ( to end trypsin activity)
- Transfer to Gimsa stain 4 mins.
- All the staining solution must be kept at 37 °C water bath during staining procedure.
- Rines in D-water
- Leave to dry completely
- Examine under microscope


# Analysis:

- ❑ Count chromosomes in 10-15 metaphases
- ❑ Count 30 if mosaicism suspected
- ❑ Detailed analysis of 3-5 metaphase

# G banding :



Normal Male Karyotype



Normal Female Karyotype



# Q banding

- ❑ Used for Y chromosome abnormalities or mosaicism
- ❑ Similar to G band ( but It can detect polymorphism)
- ❑ Examined under fluorescent microscope

# R banding:

- ❑ Used to identify X chromosome abnormality
- ❑ Heat chromosomes before treating with Gimsa
- ❑ Light and dark bands are reversed

# C banding

To identify centromere/ heterochromatin

## Heterochromatic region:

- contains highly repetitive DNA sequence
- highly condense chromatin fibres
- genetically inactive (structural elements)

## •Treated chromosomes:

- Acid
- Alkaline
- Then G band

# Karyotyping:



Organization of the chromosomes of an individual, lined up using computer image processing according to:

Size, the largest to smallest


Location of centromere

Banding pattern


# Karyotyping:



# Karyotyping:



Processed image




Karyotype

# Banding pattern:

A 'band' is defined as that part of a chromosome which is clearly distinguishable from its adjacent segments by appearing **darker or lighter** with one or more banding techniques.

# Chromosome condensation:

Various stages in condensation of DNA and Chromatin to form a metaphase chromosome






# High resolution banding:

- ❑ prometaphase and prophase chromosomes
- ❑ treatment of the cultures with chemical agent to produce cell synchrony such as (Methotrexate).
- ❑ then adding a release agent to prevent contraction such as (Thymidine).

# High resolution banding:



**Resolution G-banding:** G-banded chromosome 22 at various levels from 400-850 band level. Illustrate how a prophase and prometaphase bands coalesce to form metaphase bands.