


Blood Biochemistry BCH 471[Practical]


**Lab (4) Hemolyzing Agents & Detection of Blood**



# Blood Hemolysis

---

- **Hemolysis** (from the Greek **Hemo**: meaning blood, **lysis**, meaning to break open).
- It is the **breaking open** of red blood cells and the release of hemoglobin and the red cell contents into the surrounding fluid (plasma).
- Hemolysis may occur *in vivo* or *in vitro*.




# Hemolysis in Vivo

- **Conditions that can cause hemolysis include:**

1. Immune reactions
2. Infections
3. Medications
4. Toxins and poisons

- Because the concentration of potassium inside red blood cells is much higher than in the plasma and so elevated potassium is usually found in biochemistry tests of hemolysed blood.



# Hemolysis in Vitro

---

1. Improper technique during collection (e.g. incorrect needle size, excessive suction)
2. pH imbalance (addition acid or base)
3. Placing RBCs in a hypotonic solution

**Note:** In this lab blood hemolysis will be done by using hypotonic solutions and pH imbalance.

## When Blood Hemolysis Should Be Done?

- Breaking down RBCs to release their content
- Estimation of hemoglobin
- To obtain erythrocyte free preparation of leukocyte and platelet


# Osmosis and Osmotic Pressure

## Osmosis:

- It is the diffusion of solvent molecules across a semi-permeable membrane into a region of higher solute concentration.
- Once an equilibrium is reached the flow of water stops.

**Osmotic pressure:** the pressure exerted by a solvent passing through a semi-permeable membrane in osmosis.

**Tonicity:** the concentration of a solution as compared to another solution.




# Tonicity

## Types of solutions:

### ➤ Isotonic

- A solution that has the same solutes concentration as the normal cells of the body and the blood, having equal **osmotic pressure**.
- Example of Isotonic solution is **sodium chloride 0.9%**, have the same osmotic pressure as serum and they do not affect the membranes of the RBCs.
- In hospitals, intravenous fluids are isotonic.

Solute inside the cell = Solute outside the cell



# Tonicity



---

## Types of solutions:

### ➤ Hypotonic

- In a hypotonic solution, there is a lower concentration of solute outside a cell, creating an environment with lower osmotic pressure than what is contained within the cell.
- The RBCs will burst or hemolyzed.
- Any concentration of NaCl that is **lower than 0.9%**, will be considered hypotonic for cells.

Solute outside the cell < Solute inside the cell




# Tonicity

## Types of solutions:

### ➤ Hypertonic

- In a hypertonic solution, there is a higher concentration of solute outside a cell, creating an environment with higher osmotic pressure than what is contained within the cell.
- The RBCs will be shrink.
- Any concentration of NaCl that is **higher than 0.9%**, will be considered hypertonic for cells.

Solute outside the cell > Solute inside the cell





# Practical Part

---

# Objectives

---

1. To detect the presence of hemolysis in blood sample.
2. To detect the presence of blood in a biological sample.

# Calculations

## How to Calculate the Concentration of an Isotonic Solution of a Specific Substance:

- For example you want to know the grams of NaCl that will make a 100 ml isotonic solution, knowing that the osmolarity of RBC = 0.308 Osmolar.

### First: Calculate the molarity from osmolarity equation: [1]

Osmolarity = 0.308 Osmolar

No. of dissociation particles = 2, since  $\text{NaCl} \rightarrow \text{Na}^+ + \text{Cl}^-$

$$\rightarrow M = \frac{\text{Osmolarity}}{n} = \frac{0.308}{2} = 0.154 \text{ M}$$

$$[1] \text{ Osmolarity} = M \times n$$

Where:

M = molarity

n = No. of dissociation particles

## Calculations

**Second:** Calculate the No. of moles expressed in (w/v %): [2]

To calculate in w/v %  $\rightarrow M = \text{No. of moles} / V \text{ (in L)}$

$\rightarrow \text{No. of moles} = M \times V \text{ (in L)} =$

$\rightarrow 0.154 \text{ (from step 1)} \times 0.1 \text{ (100 ml, because you want it as \%)} = 0.0154 \text{ moles}$

**Third:** Calculate weight in grams knowing that Mwt of NaCl = 58.5 g/mol: [3]

$\rightarrow \text{Wt (g)} = \text{No. of moles} \times \text{Mwt} =$


$\rightarrow 0.0154 \text{ (from step 2)} \times 58.5 = 0.9 \text{ g}$

$= 0.9 \% \rightarrow$  the concentration of NaCl that will make an isotonic solution

$$[2] \text{ Molarity} = \frac{\text{No. of moles of solute}}{\text{Volume (L)}}$$

$$[3] \text{ No. of moles} = \frac{\text{Wt(g)}}{\text{Mwt}}$$

# Experiment (1): Hemolysis Test



## Method

1. Label 6 tubes (A → F). Then, add 1ml of RBCs suspended in saline into each tube


	Tube A	Tube B	Tube C	Tube D	Tube E	Tube F
NaCl 0.45%	5 ml					
NaCl 1.2%		5 ml				
Sucrose 6%			5 ml			
NaOH 0.1M				3 drops		
HCl 0.1 M					3 drops	
Dis. Water						5 ml
NaCl 0.9%				5 ml	5 ml	

2. Wait 30 min
3. Observe wither hemolysis has taken place



Pause and Think **What type of solution is distilled water considered?**

# Results



**A** Normal, non-hemolyzed sample

**B** Sedimented after one hour

**C** Hemolyzed sample

**Note:** the hemolyzed sample is transparent, because there are no cells to scatter the light.




## Experiment (2): Detection of Blood by Benzidine Test

---


- It is often necessary to detect the presence of small quantities of blood in urine, stomach contents etc.

### Principle

- This method depend on the fact that the **heme group of hemoglobin** possesses a peroxidase-like activity which catalyzes the breakdown of hydrogen peroxide ( $\text{H}_2\text{O}_2$ )
- The oxidizing species formed in this reaction can then react with benzidine giving **blue greenish color**.



**Note:** the test is **not specific** for blood as peroxidases present in milk, potatoes and pus, as well as the ions of  $\text{Fe}^{+3}$ ,  $\text{Cu}^{+2}$  and  $\text{K}^{+1}$  will give false positive results



# Experiment (2): Detection of Blood by Benzidine Test

---

## Method

- Place 3ml of sample in a boiling water bath for 3 minutes.
- Cool it under tap water.
- Add 2 ml Benzidine+ 1 ml H<sub>2</sub>O<sub>2</sub>

## Results

- If the test is negative → blood is absent from sample.
- If the test is positive → blood is probably **not definitely** present in sample.
- For this reason these tests are often described as **“presumptive tests”** .



Positive results

