

Polymerase Chain Reaction (PCR)

By: Sahar AlSubaie

PCR:

PCR is a technique for amplifying DNA sequences in vitro.


PCR properties:

- Rapid & easy.
- Sensitive.
- Robust.
- Widespread applications.

Basic requirements for PCR reaction:

- DNA sequence of target region
- Primers
- Taq polymerase (Thermo-stable DNA polymerase)
- dNTPs
- DNA thermal cycler

Basic requirements for PCR reaction:



Standard PCR reaction:

- DNA template (target) to be amplified. 1-5uL.
- Two primers (forward and reverse).
- Master mix.

Master mix consist of :

- 1-2ul Taq polymerase. (optimum temperature at around 70C.)
- 1-2.5ul dNTPs.
- 4-7ul Buffer solution
- 1-2ul Divalent cations (Mg^{2+})
- 1-2ul Monovalent cation (K^{+})

PCR cycling steps:

- Denaturation
- Annealing
- Extension

Denaturation:

Double stranded DNA are separated to form single stranded DNA.

Temperature: 93-95C°.

Annealing:

Primers bind to their complementary sequences.

Temperature: 50-70C° (depending on the melting temperature of the expected duplex).

Extension:

DNA polymerase binds to the annealed primers and extends DNA at the 3' end of the chain.

Temperature: 72C°

PCR cycling steps:

PCR : Polymerase Chain Reaction

30 - 40 cycles of 3 steps :

Step 1 : denaturation


1 minut 94 °C



Step 2 : annealing


45 seconds 54 °C

forward and reverse primers !!!



Step 3 : extension

2 minutes 72 °C
only dNTP's



(Andy Vierstraete 1999)

Standard PCR method:

- ❑ Prepare mix containing primers, dNTP's, buffer, Taq polymerase and water sufficient for all reaction tubes.
- ❑ Aliquot appropriate volume to each tube.
- ❑ Add DNA to each tube using a new tip for each sample.

Standard PCR method:

- Load tubes on PCR machine.
- Cover the plate.

N.B If PCR machine does not have heated lid, 1 drop of mineral oil should be added to each tube before loading then start the program.

Properties of Polymerase:

- ❑ Taq polymerase originally isolated from *Thermus aquaticus*.
- ❑ Heat stable (half life of about 30 min at 95C°).
- ❑ Taq DNA polymerase has no proof-reading function in 3' to 5' direction.
- ❑ Primer extension occurs at up to 100 bases/sec.
- ❑ Plateau is eventually reached.

Trouble Shooting:

Problem	Possible reasons
No product	Primers annealing?
Product of incorrect size	Primers annealing elsewhere in genome?
Several Products formed	Contamination? Several annealing sites?

PCR Rxn inhibitors:

- ❑ proteinase K: digest polymerase
- ❑ phenol: denature the polymerase.
- ❑ EDTA: chelates Mg^{+2}

Advantages of PCR:

- Uses less patient DNA.
- Result obtained more quickly about 3hr for PCR.
- No need to use radioactive material for PCR.
- Precise in determining sizes of alleles.
- Can be used to detect point mutations.

Disadvantages of PCR:

- ❑ Target DNA sequence must be known.
- ❑ Errors made by Taq polymerase.
- ❑ Size limitations particularly in GC rich triplet repeats.

PCR Based Technologies:

- Multiplex PCR
- QF-PCR
- RT-PCR