

# ORAL GLUCOSE TOLERANCE TEST (GTT)

**- Objectives:**

- Use OGTT in diagnosis of diabetes mellitus.

## - INTRODUCTION

- Serial measurement of plasma glucose before and after glucose is given orally should provide a standard method to evaluate individuals and establish values for normal and disease states.
- There are two type of glucose tolerance test (Oral and IV).
- The most common glucose tolerance test is **the oral glucose tolerance test (OGTT)**.
- The test reveals how quickly glucose is metabolized from the bloodstream for use by cells as energy source.

- There are a number of factors that may affect glucose tolerance and that should be controlled or eliminated before such a test is performed :

cigarettes

Anxiety

activity

Coffee

amount of glucose ingested

carbohydrate intake

Time of pervious food intake

corticosteroids

Age

inactivity

weight

Some medicines

## - How the test is performed

- When an oral glucose tolerance test is ordered, the following conditions should be met:

- (1) Omit medications known to affect glucose tolerance.
- (2) Perform the test in the morning after 3 days of unrestricted diet and activity.
- (3) Perform the test after a 10-16 hours fast (**better 12 hour**).

- **Oral dose** : For adults, the recommended load is 75 g and for children, 1.75 g/kg,

- Plasma glucose should be measured **fasting** then every 30 min for 2h **after** an oral glucose load

- **Note**: the time of collection is different, it is depend on the situation.

## - Why might I need to have the test?

- Generally most healthcare providers recommend that all pregnant women be screened for gestational diabetes.
- Experts recommend this test to pregnant women who are between 24 and 28 weeks of pregnancy .
- This test is also recommended for anyone suspected of developing adult diabetes.

## **- How the side effect during the test?**

- Some people feel sweaty, light-headed, or may even feel short of breath or faint after drinking the glucose.
- However, serious side effects of this test are very uncommon.

## Normal and abnormal results :

**Normal :** Fasting: 60 -128 mg/dL  
1 hour: less than 200 mg/dL  
2 hours: less than 140 mg/dL

## **Abnormal :**

-Higher-than-normal levels of glucose may mean you have **prediabetes** ,  
**diabetes (type 2)**, or **gestational diabetes**.

-Between 140 - 200 mg/dL is called *impaired glucose tolerance*. And this  
Called "prediabetes." It means you are at increased risk for developing  
diabetes.


- A glucose level of 200 mg/dL or higher is a sign of diabetes (in adult  
individual) or gestational diabetes (in pregnant woman).

-**However**, high glucose levels may be related to another medical problem.



## - PRINCIPLE:

- Several reagents can be used to assay reducing sugars such as 3, 5 dinitrosalicylic acid in one of the compounds.
- In alkaline solution it is reduced to 3-amino-5- nitro salicylic acid, which is orange-red.
- Absorbance is determined at 540 nm.



## - GLUCOSE ESTIMATED BY O-TOLUIDINE METHOD:

|                            | Plasma | Standard | dH2O | DNS reagent |
|----------------------------|--------|----------|------|-------------|
| Test (a1) (Fasting plasma) | 0.1    | -        | -    | 2 ml        |
| Test (a2) (Fasting plasma) | 0.1    | -        | -    | 2 ml        |
| Test (b1) Tow- hour        | 0.1    | -        | -    | 2 ml        |
| Test (b2) Tow- hour        | 0.1    | -        | -    | 2 ml        |
| Standard (1)               | -      | 0.1      | -    | 2 ml        |
| Standard (2)               | -      | 0.1      | -    | 2 ml        |
| Blank                      | -      | -        | 0.1  | 2 ml        |

Mix the contents of each tube and cover each tube by Aluminum foil

Boiling water bath for 5 minutes

cool the tubes for 1-3 min

Read absorbance at 540 nm

## - RESULT:

| <b>Tubes</b> | <b>Absorbance at 540 nm</b> |
|--------------|-----------------------------|
| Test (a1)    |                             |
| Test (a2)    |                             |
| Test (b1)    |                             |
| Test (b2)    |                             |
| Standard (1) |                             |
| Standard (2) |                             |

## - CALCULATIONS:

- Conc. Of Std. = 0.1 g/dl.

- Sample A = Fasting plasma glucose

- Sample B = Two hour plasma glucose

- Amount of glucose in plasma =  $\frac{\text{Means Ab Test}}{\text{Means Ab Std.}} \times \text{conc. Of Std} = Z \text{ g/dl}$

-  $Z \text{ g/dl} \times 1000 = \underline{Y \text{ mg/dl}}$

- Calculate the glucose in fasting glucose plasma and in two hours plasma glucose ..

- Then discuss your results ..