




# **Triglyceride determination**



## Introduction:

- Triglycerides are esters of fatty acids and are hydrolyzed to glycerol and free fatty acids (by lipase).
- Triglyceride determinations when performed in conjunction with other lipid assays are useful in the diagnosis of **primary and secondary hyperlipoproteinemia**.
- They are also of interest in following the course of diabetes mellitus, nephrosis, biliary obstruction, and various metabolic abnormalities due to endocrine disturbances.







- **Hyperlipoproteinemia:** abnormally elevated of fat in blood ( disorder in lipid metabolism).

- Standard methods for the measurement of triglyceride concentrations involved either **enzymatic** or **alkaline hydrolysis** to liberate glycerol.

## - Principle:

The enzymatic reaction sequence employed in the assay of Triglycerides is as follows:




- The present procedure involves hydrolysis of triglycerides by lipase .

- **The glycerol concentration** is then determined by enzymatic assay coupled with **Trinder reaction** that terminates in the formation of **a quinoneimine dye** .

- **The amount of the dye formed**, determined by its absorption at 500 nm, is directly proportional to the concentration of triglycerides in the samples.

## - Specimen collection and storage:

1. Fresh, non-hemolyzed serum from fasting patients is recommended.
  2. Triglycerides in serum appears stable for three days when stored at 2-8 °C.
  3. Prolonged storage of the samples at room temperature is not recommended since other glycerol containing compounds may hydrolyze, releasing free glycerol with an apparent increase in total triglycerides content.
- 

## - Method:

- By Triglyceride reagent kit.

-Follow the table:

	Blank	Standard	Test
Reconstituted Reagent	1 ml	1 ml	1 ml
<b>Pre-worm at 37°C for 2 min and add:</b>			
Distilled water	0.01 ml (10 µl)		
Standard	---	0.01 ml (10 µl)	
Sample	---	---	0.01 ml (10 µl)
<b>Mix and incubate at 37°C for 5 min</b> ↓ <b>Read the absorbance of standard and sample at 500 nm against blank</b>			

## **-Calculation:**

$$\text{Conc .of TG} = \frac{\text{Ab Test}}{\text{Ab Std.}} \times \text{conc. of Std. ( 200 mg/d)}$$

**- Normal range:** 10 -190 mg/dl




# **HDL-Cholesterol determination**






## - Introduction:

- Cholesterol is a fatty substance found in blood, bile and brain tissue.
- It serves as a precursor to **bile acids, steroids and vitamin D**.
- **In the plasma, cholesterol is transported by three lipoproteins:** high density lipoprotein (HDL-Cholesterol), low density lipoprotein (LDL-Cholesterol), and very low density lipoprotein (VLDL- Cholesterol).
- The concentration of **total cholesterol** in serum has been associated with metabolic, infectious and coronary heart diseases.



- The concentration of HDL-cholesterol in serum has important in diagnosis of the how the level of **risk to get coronary heart diseases**.
  - More HDL-Chol. That indicate low risk to get coronary heart disease.
  - Castelli and co-workers have indicated that an **inverse relationship** exists between serum HDL-Cholesterol and the risk of coronary heart disease.
  - The measurement of **HDL Cholesterol and triglyceride** provides valuable information for the prediction of coronary heart disease and for lipoprotein phenotyping.
- 



## - Specimen collection:

1. Specimen should be serum and free from hemolysis.
2. Patient should be fasting for 12-14 hours.

## - Principle:

### - HDL cholesterol determination

- It is direct method without specimen pre treatment.



- A color reaction which is proportional to HDL-Cholesterol concentration. **The absorbance is measured at 600 nm.**

**LDL**= Low density lipoprotein

**POD**= Peroxidase

**VLDL**= Very low density lipoprotein

**CE**= Cholesterol Esterase

**CO**= Cholesterol Oxidase

**AAO**= Ascorbate Oxalase

**4-AAP**= 4- Aminoantipyrine

**DSBmT**= N,N-bis (4-sulphobutyl)-m-toluidine-disodium

## Method:

### - HDL Cholesterol:

- Follow the Table:

	Blank	Calibrator	Assay
Reagent R1	300 $\mu$ l	300 $\mu$ l	300 $\mu$ l
Calibrator		3 $\mu$ l	3 $\mu$ l
<b>Mix vigorously, let stand for 5 min at 37°C . Read absorbance A1 at 600 nm against blank.</b>			
Add Reagent R2	100 $\mu$ l	100 $\mu$ l	100 $\mu$ l
<b>Mix vigorously, let stand for 5 min at 37°C . Read absorbance A2 at 600 nm against blank.</b>			

## - Calculation:

\* Determine the HDL Cholesterol conc.

$$\Delta \text{ Abs.} = (A2 - 0.75 A1)$$

$$\text{Conc. Of HDL} = \frac{\Delta \text{ Ab Assay}}{\Delta \text{ Ab Calibrator}} \times \text{conc. of calibrator (50 mg/dl)} = \text{mg/dl}$$

## - Normal value of:

### - HDL-Cholesterol :

- Low level (risk factor) < 40 mg/dl

- High HDL (protector factor)  $\geq$  60 mg/dl