

Alternative Modes of Mechanical Ventilation

Dr. Zohair Alaseri, MD
FRCPC, Emergency Medicine
FRCPC, Critical Care Medicine
Intensivist and Emergency Medicine Consultant
Chairman, Department of Emergency Medicine
King Saud University Hospitals, Riyadh, KSA

Alternative Modes of Mechanical Ventilation

Dual-control ventilation modes

Advantages

volume-control
(guaranteed minute
ventilation)

+

pressure-control ventilation
(limited peak airway
pressure).

- Increase the safety and comfort of mechanical ventilation.
- No RCT indicate improved patient outcomes (including mortality).

Alternative Modes of Mechanical Ventilation


Pressure-Regulated Volume-Control (PRVC)

- This mode is under the dual control of pressure and volume.
- The physician presets a desired TV, and the ventilator delivers a pressure-limited (controlled) breath until that preset TV is achieved.
- The breath is essentially like a conventional pressure-controlled ventilation breath, but the ventilator can guarantee a predetermined minute ventilation.

Alternative Modes of Mechanical Ventilation

Pressure-Regulated Volume-Control (PRVC)

- Breath to breath, the inspiratory pressure is automatically adjusted down or up to deliver a preset TV.
- If the delivered volume is too low, it increases the inspiratory pressure on the next breath.
- If it is too high, it decreases the inspiratory pressure.
- This adjustment gives the patient the lowest peak inspiratory pressure needed to achieve a preset TV.
- Advantage.....deliver minimum minute ventilation at the lowest peak airway pressures possible



The pressure, volume, and flow to time waveforms for pressure-regulated volume-controlled ventilation


Alternative Modes of Mechanical Ventilation

Dual-control breath-to-breath, pressure-limited, flow-cycled ventilation

- volume-support ventilation (VSV) or variable-pressure-support.
- Combination of PSV and volume-control ventilation.
- Like PSV, the patient triggers every breath, controlling his or her own respiratory frequency and inspiratory time.
- This mode delivers a breath exactly like conventional PSV, but the machine can guarantee minute ventilation.
- The pressure support is automatically adjusted up or down to deliver a preset TV.

Alternative Modes of Mechanical Ventilation

Volume-support ventilation (VSV) or Variable-pressure-support

- it is flow cycled, which means that the patient determines the respiratory rate and inspiratory time.
- The mode cannot be used in a patient who lacks spontaneous breathing effort.
-  Volume support has also been marketed as a self-weaning mode.

Alternative Modes of Mechanical Ventilation

Volume-support ventilation (VSV) or Variable-pressure-support

Potential Problems

- If the patient's metabolic demand increases, raising the tidal volume, the pressure support decreases to provide less ventilatory support when the patient needs it most result in hypoxemia.
- TV must be correctly set to the patient's metabolic needs.
- If the tidal volume is set too high, weaning is delayed. If it is set too low, the work of breathing increase

Alternative Modes of Mechanical Ventilation

3--Automode and variable support or variable-pressure control

- This mode is basically the combination of the 2 modes described above.
- If the patient has no spontaneous breaths, the ventilator is set up in the PRVC mode.
- When the patient takes 2 consecutive breaths, the mode is switched to **Volume-support ventilation VSV**.

Alternative Modes of Mechanical Ventilation

3--Automode and variable support or variable-pressure control

- Designed for automatic weaning from pressure control to pressure support depending on the patient's effort.
- no randomized trials have been conducted to show this type of weaning is more effective than conventional weaning.

Alternative Modes of Mechanical Ventilation

Dual control within a breath

- This mode has been called volume-assured pressure support or pressure augmentation.
- This mode can switch from pressure control to volume control within a single specific breath cycle. After a breath is triggered, rapid and variable flow creates pressure to reach the set level of pressure support.

Alternative Modes of Mechanical Ventilation

Dual control within a breath

- TV is monitored.
- If it equals the minimum set TV, the patient receives a typical pressure-supported breath
- if the TV is less than the set one, the ventilator switches to a volume-controlled breath with constant flow rate until the set tidal volume is reached.

Alternative Modes of Mechanical Ventilation

Dual control within a breath

One study compared volume-assured pressure support with simple assist-control volume support and showed a 50% reduction in the work of breathing, lowered airway resistance, and lowered intrinsic PEEP.

Alternative Modes of Mechanical Ventilation

Automatic tube compensation

- Specifically for weaning
- Designed to overcome the resistance of the endotracheal tube by means of continuous calculations.
- These calculations allow the ventilator to supply the appropriate pressure needed to overcome this resistance
- no studies so far

Alternative Modes of Mechanical Ventilation

Proportional assist ventilation


- Decrease the work of breathing.
- The mode adjusts airway pressure in proportion to the patient's effort.
- lets the patient determine the inspired volume and the flow rate.
- The support given is a proportion of the patient's effort and is normally set at 80%.

Alternative Modes of Mechanical Ventilation

Proportional assist ventilation

- This support is always changing according to patient's effort and lung dynamics.
- The patient's work of breathing remains constant regardless of his or her changing effort or demand
- This mode can be used only in patients with spontaneous respiratory efforts.
- Not approved by FDA

Proportional assist ventilation PAV




Alternative Modes of Mechanical Ventilation

Airway pressure–release ventilation

- Bilevel, or biphasic, ventilation
- new mode.
- The ventilator is set at 2 pressures
 high CPAP & low CPAP both levels are time
 cycled.

APRV



Alternative Modes of Mechanical Ventilation

Airway pressure–release ventilation

- The high pressure is maintained for most of the time, while the low pressure is maintained for short intervals of usually less than 1 second to allow exhalation and gas exchange to occur.
- The patient can breathe spontaneously during high or low pressure

Alternative Modes of Mechanical Ventilation

Airway pressure–release ventilation

- Has benefit of alveolar recruitment.
- Its disadvantage is that the tidal volume is variable.
- The clinician must be constantly aware of the patient's minute ventilation to prevent severe hypercapnia or hypocapnia.