

BCH302
Second semester 1432-1433



Student name:

Student Number:

Question(1):

Put true "T" or False "F" and correct the false one:

1-The lead sulfite test is specific for Methionine []

.....

2-Each protein can precipitate at different salt conc. depend on their Molecular weight. []

.....

3-A homo-polymeric saccharide is a polysaccharide consist of mixture of different monosaccharide. []

.....

4- In Iodine Test the Iodine forms a large blue complex with Alfa- amylase helix, []

.....

5- when calcium chloride added to the saponification solution the soap get out of solution and separated in the form of insoluble floats above the surface []

.....

6- Iodine reacts with unsaturated compounds, The double bond is changed to single bond and iodine is added to both side of the bond []

.....

7- The Millon's test is specific for only amino acid contain hydroxyl group attached to benzene ring []

.....

8-The denaturation of protein lead to loss the primary structure of protein []

.....