

# Quiz

Name:-

Student number:-

- 1- In the SDS-PAGE gel there are two different layers, the upper layer called.....
- 2- In the western-blot experiment the primary antibody must be specific to the.....
- 3- In the western-blot a stain can be optionally used called Ponceau S, which gives ..... colour when there is a protein sample, and this stain is removable by washing by water.
- 4- In SDS-PAGE experiment a tracking dye called .....is used to track the sample in the gel.
- 5- There are several stains can be used to stain the samples in the SDS-PAGE experiment, such as.....