

إجابة نموذجية لأحد الاسئلة القصيرة
(Quiz)
Quiz

Name:-

Student number:-

1- In the SDS-PAGE gel there are two different layers, the upper layer called.....*Stacking gel*

2- In the western-blot experiment the primary antibody must be specific to the.....*targeted protein.*

3- In the western-blot a stain can be optionally used called Ponceau S, which gives*Red*..... colour when there is a protein sample, and this stain is removable by washing by water.

4- In SDS-PAGE experiment a tracking dye called*Bromophenol Blue*.....[↓]is used to track the sample in the gel.

5- There are several stains can be used to stain the samples in the SDS-PAGE experiment, such as.....*Coomassie brilliant blue.*